



## **TECHNICAL SERVICE BULLETIN # 10**

**Product:** Fabius GS

**Date:** March 29, 2004

**Re:** Software Release – 2.06

**Reference Doc:** Product Bulletin #164, SP00297

**Reason/Rational:** Enhancement – See Product Bulletin #164 for details

**Solution/Procedure:** Perform Software 2.06 upgrade on all Fabius GS units currently containing 2.n version. Complete the Return Form (attached) and return to DMI/DS within 7 days of account completion.

**Devices affected:** All Fabius GS containing 2.n **only**.

**Timing/When:** At customer request, next service or Periodic Manufacturer's Service

**Additional Info:** Fabius GS units containing software 2.n shall be updated to 2.06 by ordering Software Kit – 4118735

Note Do not install software 2.06 on Fabius GS units containing 1.3n. See TSB #9 for details.

**Cost:** Software upgrades to 2.06 are available to Fabius GS customers at no charge.

**Ordering Info/Parts:** Contact your DMI Customer Service Representative to order the software.

**Distribution:** DrägerService N.A. & all Authorized Service Organizations.

## Return Form - Software 2.06 Upgrade

Distributor:	Customer:
Country:	Address:
	City:

Machine S/N

[illegible][illegible][illegible]

Date of Update \_\_\_\_\_  
Signature (TSR) \_\_\_\_\_  
Name (Print) \_\_\_\_\_